DE/EN

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**INSTRUCTIONS  
FOR USE**

# AssistX Call

**C.S.S.**  
MICRO SYSTEMS

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
## 1. Introduction

Thank you for choosing a high quality CSS product. Please read these instructions for use carefully and take care of these instructions when using the device. Please keep the instructions for use at hand when you use the device.

## 2. Intended use

**AssistX Call** is a sensor adapter for nurse call systems and enables calls to be triggered using special sensors that cannot be connected directly to the call system. For example, with the **AssistX Call** proximity sensor **AssistX ProX**, acoustic sensor **AS one**, sipp and puff sensor Hy-Pneumax, infrared receiver **AssistX IR one** or various buttons or sensors can be used to call for help.

**AssistX Call** is tailored to the call system and delivered ready to use with an appropriate connection cable.

---

### 3. Explanation of symbols

Safety instructions!

Important instructions!

---

## 4. Safety instructions

### Handling

Please handle the device carefully. The device contains sensitive electronic components. It is not allowed to do any modification. The device is not life-sustaining or life-supporting. The device consists of several components. Small items may be detached from the device, which may result in suffocation or other danger to children. Small children or people with intellectual disabilities should not reach or use the device without supervision.

Do not use the device in damp, wet, very cold or very hot environments (see temperature information in technical data).

Never replace the built-in batteries by yourself! The battery can be damaged, which can lead to overheating and can lead to injuries. The integrated rechargeable lithium polymer battery may only be replaced by authorized personal.

### Mounting

The device must be mounted in a way that it cannot fall down. Only recommended mounting components from the manufacturer are to be used. The user bears the sole risk of mounting.

---

The manufacturer is not liable for personal injury or property damage caused by the device or its components falling down.

### **Repair**

Under no circumstances you should open the device and repair it by yourself. Disassembling the device can permanently damage the device and cause injury. Contact your reseller or the manufacturer if the device is damaged or does not work properly. If the device or its accessories are mechanically damaged, they may no longer be used.

### **Interaction with medical devices**

The device and its components can emit electromagnetic radiation and have integrated magnets. The device must therefore not be placed on a person or brought closer than 15 cm to people using other active medical devices. The device must also not be placed directly on active medical devices. Please pay special attention to people with implanted pacemakers. If necessary, please contact the device manufacturer of the other active medical device.

---

## Attention! Magnets

The mounting plate contains a strong neodym magnet. Magnets can affect cardiac pacemakers and implanted defibrillators.

- The pacemaker could be switched into testmode and cause discomfort
- A defibrillator could stop working
- Keep enough distance to magnets if you use such devices
- Advice other users of such device against getting too close to magnets
- Keep magnets away from all devices and objects, which can be damaged by strong magnetic fields

**Recommended safety distance is 6 inch (15 cm).**

Magnets cause a far-reaching magnetic field. They could damage TV sets, laptops, PC harddrives, creditcards, ATM cards, data storage, mechanical watches, hearing aids and speakers.

---

Vehicle keys, USB sticks, cameras, cell phones and smart phones contain non-magnetic storage media. Therefore, static magnetic fields near those devices cannot delete data. There are countless covers with magnetic closures or magnetic holders on the market, which also speaks against the notion that magnets can damage these devices. However, it cannot be ruled out that very strong magnetic fields might magnetise and maybe damage mechanical parts or built-in speakers. When in doubt, keep these devices away from strong magnets.

## 5. Cleaning and disinfection

Before cleaning the device, switch it off and disconnect all cables. The device should be cleaned if required with a disinfectant and cleanser intended for use in hospitals. Cleaning should only be done with a wipe. Wipe all surfaces and let them dry or follow the wiping instructions.

**Do not use sprays** that may allow moisture to enter the device. Under no circumstances may any component of the device be subjected to steam or wet sterilization.

---

## 6. Scope of delivery

- **AssistX Call**
- connection cable
- optional 5V USB power supply (depending on the call system)
- Magnetic wall mount with screws and velcro
- Instructions for Use

## 7. Device Buttons, Ports, Display

---

- 1 2x 3,5 mm jack plugs (**AssistX Port**)
- 2 Button with LED
- 3 Status LED
- 4 USB connector for optional power supply
- 5 for cable to nurse call

### Battery charge status (if present)

Briefly press the button.

	red:	< 25%
	orange:	> 25%
	light green:	> 50%
	green:	> 75%

### Status LED / Button LED

<i>Function</i>	<i>LED</i>	<i>Button LED</i>	<i>Signal</i>
<b><i>New sensor detected</i></b>	 1 sec		
<b><i>Sensor intentionally unplugged</i></b> Press the button to confirm intentionally unplugging.			

<p><b>Sensor unplugged</b> Automatic Call</p>			 when set
<p><b>Defective cable detected</b> Automatic Call</p>			 when set
			
<p><b>Sensor with normally open contact plugged in</b> Operate sensor to check function. (LED turns off after operation)</p>			
<p><b>Sensor with normally open contact not operated for more than 24h</b> Operate sensor to check function. (LED turns off after operation)</p>			
<p><b>Call suppression during checking sensor function</b> Test mode is active whenever an unsafe</p>			

<p>sensor has been plugged in or a function test has been requested by pressing the button.</p> <p>Test mode finishes if an operation of every unsafe sensor has been detected. After 20 sec test mode finishes anyway.</p>			
<p><b>No external power supply Device powered by battery</b> Automatic call once if the external power supply fails.</p>			
<p><b>Low battery</b></p>			
<p><b>External powersupply connected (again)</b></p>			

## 8. Use of the device

**AssistX Call** is supplied with power via the call system or an external 5V USB power supply unit. In both cases, depending on the call system, the device also works in the event of a power failure in the network. When powered by an ④ external power supply unit, the built-in lithium polymer battery guarantees full functionality for a few hours in the event of a power failure.

Depending on the call system, the call line can be interrupted or closed for an adjustable time. The switching contacts are connected to the call system using the cable that matches the call system.

### **Standard configuration:**

In the standard configuration, **AssistX Call** switches the call contact for 1.5 seconds when the entrance sensor has been activated.

- **Connect Sensors**

Connect one or two sensors to ① sockets S1 and S2 as required.

---

- **Establish connection**

Establish the connection to the call system at the ⑤ RJ45 modular socket with the enclosed 8-pin cable.

### **Switching ON and OFF:**

If the device is not required for a longer period of time, the electronics can be switched off. The device switches on automatically when it is connected via the USB socket or an external power supply is connected through the call system.

To switch on, press the ② button for 2 sec. Release the ② button after the tone signal, the ③ LED lights up **green**. To switch off, press the ② button for 3 sec. Release after the tone signal. After the flashing, the ③ LED turns off, the ② button LED turns **red**.

After the device has been switched on, the sensor inputs are checked. The device automatically adjusts itself to the connected sensors. If single sensors with a closing contact are plugged in, their functionality must be checked. The ③ LED in the button flashes until each sensor with a closing contact has been actuated once. The call release is suppressed for 20 sec during the test period.

---

## **Feedback:**

The caller gets informed by the device when the call release has been recognized and is forwarded. Two different „Call recognition“ situations are treated. On the one hand, the dispatch of the call by the **AssistX Call**, on the other hand, the feedback from the call system that the call has been received.

- ***Call sent***

The sending of a call is signalled optically and acoustically. The signals inform, that **AssistX Call** has recognized the sensor action and switched the call line. The signal is shown regardless of the call system having recognized the call or not.

- ***Confirmation indicator***

Call buttons must have a so-called confirmation lamp, which is controlled by the call system to signal the call triggering to the caller. The confirmation light is **red**. **AssistX Call** receives the confirmation light signal from the call system and uses it to switch on the **red** **③** LED. This signal is controlled almost directly by the call system and indicates that the call has been received.

---

## Control call:

The **AssistX Call** monitors safety-relevant parameters. If proper function is no longer guaranteed, the device triggers a call. Be monitored:

- **Power supply for the control logic**  
A simple call is triggered
  - **ensor error ejection**  
Optical/acoustical call
  - **Sensor error cable break (for sensors with NC)**  
Optical/acoustical call
  - **Function check (for sensors with NO)**  
Optical call **red** ③ LED after 24 hours without actuation,  
Actuation of the sensors required
  - **Sensor test as required (for sensors with NC)**  
Press button (settings)
  - **Sensor test on request (for sensors with NO)**  
Press button (settings)
-

- **Sensor test after plugging in (for sensors with NO)**

Flashing ② button LED, goes off when the sensors are actuated; A time window of 30s is provided for checking the sensor function. This is indicated by the **orange** ③ LED. During this time, the call triggering is suppressed if an unsafe sensor is activated. Call suppression has no effect on safe sensors.

## 9. Settings

No settings can be made. Only the following functional tests:

<i>Function</i>	<i>LED</i>	<i>Button LED</i>	<i>Signal</i>
<b>Functional test</b> <b>Simple sensors with NO</b> Press the button for 1 sec. Release after flashing for <b>test mode</b> .			
<b>Test mode</b> <b>Simple sensors with NO</b> 20 sec active. Operate sensor to check function. The ③ LED turns out when each sensor with a closing contact has been actuated once. Checked for 24h.			

<p><b>No Simple sensors with NO connected</b> No manual sensor test necessary. ② Button LED goes out after test mode.</p>			
<p><b>Stand by mode</b> Press the button for 3 sec. Release after the signal tone. The ③ LED turns off, the ② button LED remains <b>red</b>.</p>			
<p><b>Turn on the device</b> Automatically or by pressing the button for 2 sec. Release after the signal tone.</p>			

## 10. Maintenance / Review

The following must be checked on the device every three months:

- Functional test
- Connectors for damage and soiling
- Damage to the housing
- If a battery is present: exchange every 3 years

## 11. Support

If you need any support or for technical questions about the product, please contact the reseller where you purchased this product. If necessary, e. g. the reseller is not available or cannot help, send us an email to [support@csslabs.de](mailto:support@csslabs.de). Please be sure to include a description of the problem and the serial number of the device.

---

## 12. Technical Data

**Dimensions:** 74 mm x 74 mm x 27 mm

**Weight:** 100 g

**IP Rating:** IP 54

**Power supply:** via call system 6.5- 24V DC, if required: external power supply via mini-USB, 5V DC, max.500mA with bridging of power failure with LiPo battery min. 24h

**Power consumption:** max. 200 mA at 5V, max. 100 mA at 12V, max. 50 mA at 24V, with ext. Power supply/battery backup max. 500mA

**connection:**

Switch connectors: 2x 3.5 mm jack plugs, **AssistX Port**

Call system: for special cable to the nurse call system

---


## 13. Conformity/Standards

 This product is CE compliant.

**AssistX Call** meets the requirements of DIN VDE0834.

## 14. Environmental protection and disposal

All electrical and electronic devices must not be disposed with general household waste, but must be disposed at designated disposal points. This product is subject to the European Directive 2002/96/EC and must therefore be properly disposed at the WEEE collection points. werden.

WEEE DE94635635

---

**THE PASSIONATE ASSISTIVE TECHNOLOGY MANUFACTURER**

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© 2024 CSS MicroSystems GmbH