DE/EN (EU)

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**INSTRUCTIONS
FOR USE**

BIF8 plus

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
WEEE DE94635635

Table of Contents

1. Introduction	15
2. Intended use	15
3. Explanation of symbols	15
4. Safety instructions	16
5. Cleaning and disinfection	17
6. Scope of delivery	17
7. Device Buttons, Ports, Display	18
8. Use of the device	19
9. Settings (only by trained specialists)	21
10. Maintenance / Review	23
11. Support	23
12. Technical specifications	24
13. Conformity/Standards	24
14. Environmental protection and disposal	24

1. Introduction

Thank you for choosing a high quality CSS product. Please read these instructions for use carefully and take care of these instructions when using the device. Please keep the instructions for use at hand when you use the device.

2. Intended use

The Bed Interface **BIF8plus** enables the wireless operation of an electrically adjustable bed with the help of an infrared remote control or an environment control device. The **BIF8plus** is plugged into the motor control of the bed instead of or with the wired hand control. **BIF8plus** provides up to ten different switching functions for this purpose.

3. Explanation of symbols

 Safety instructions!

 Important instructions!

4. Safety instructions

Handling

Please handle the device carefully. The device contains sensitive electronic components. It is not allowed to do any modification.

The device is not life-sustaining or life-supporting. The device consists of several components. Small items may be detached from the device, which may result in suffocation or other danger to children. Small children or people with intellectual disabilities should not reach or use the device without supervision.

Do not use the device in damp, wet, very cold or very hot environments (see temperature information in technical data).

Mounting

The device must be mounted in a way that it cannot fall down. Only recommended mounting components from the manufacturer are to be used. The user bears the sole risk of mounting. The manufacturer is not liable for personal injury or property damage caused by the device or its components falling down.

Repair

Under no circumstances you should open the device and repair it by yourself. Disassembling the device can permanently damage the device and cause injury. Contact your reseller or the manufacturer if the device is damaged or does not work properly. If the device or its accessories are mechanically damaged, they may no longer be used.

5. Cleaning and disinfection

Before cleaning the device disconnect all cables. The device should be cleaned daily with a disinfectant and cleanser intended for use in hospitals. Cleaning should only be done with a wipe. Wipe all surfaces and let them dry or follow the wiping instructions.

Do not use sprays that may allow moisture to enter the device. Under no circumstances may any component of the device be subjected to steam or wet sterilization.

6. Scope of delivery

- #20290 **BIF8plus** bed interface
- IR remote control
- Instructions for use

Please note: Scope of delivery may vary depending on bed type!

7. Device Buttons, Ports, Display

1 Case **BIF8plus**

2 Signal light

3 IR remote control for **BIF8plus**
with function keys 1-10

2 **Signal light**

- Green: Ready to use, stand by mode
 - Off: Off / Stays off until new signal arrives
 - Orange: Receiving of valid infrared command
 - Red: Error, functions blocked
 - Flashing red: Receiving command, error: no longer separating relay contact, no execution
-

8. Use of the device

Using the infrared remote control or the CSS setup software for PC, a **trained technician** can group any number of the ten functions for switching. Combined functions can be assigned to each function key on the remote control.

Installation:

Connect your **BIF8plus** to power. If there are no other **BIF8plus** devices in the room, you can now control the bed with your environment control device. If necessary, adjust the control further to suit your needs. Once everything is set up, the **BIF8plus** can be placed at your bed or close to the bed.

Standard control:

The bed interface comes with a standard IR remote control, which has the usual functions of a bed control in the first three sections. The last two sections can be assigned as individual commands **by a trained technician**.

Safety

BIF8plus has several security features.

Automatic switching off of a permanent function

A permanent function by continuously pressing a **3** function key is interrupted after a predefined time. The technician can set the predefined time before automatic switch-off from 1 to 25 sec.

Multiple activation of the same functions

In order to ensure safe operation via infrared, a mandatory function pause of 2 sec is inserted between the repeated activation of the same functions before the automatic switch-off is reached (with newer IR remote controls this is implemented via a toggle bit). When the automatic switch-off is reached, however, the same function can be activated again immediately. If a new **3** function key is pressed, the selected function (function group) is executed immediately.

Different device addresses

If several beds are used with a **BIF8plus** in a room, it must be ensured that each bed only reacts to the associated ③ remote control. For this reason, the receiving address can be set **by a trained technician**. One of four possible device addresses can be assigned to each **BIF8plus**, so that up to four **BIF8plus**-controlled beds can be in one room.

Monitoring switch contacts

BIF8plus recognizes when one of the switching contacts is no longer electrically separated due to a defect. In this case, the activation of the switching contacts is deactivated. This faulty condition is indicated by the ② red signal light on the front of the case.

System self-monitoring

Safety precautions have been taken in the **BIF8plus** to prevent the control unit from malfunctioning. In the event of a failure, a watchdog module deactivates the activation of the switching contacts.

High protection class

The electronics are protected against water with a housing of protection class IP65. The power supply is supplied via a short-circuit proof power supply according to safety standard IEC 60601-1.

9. Settings (only by trained specialists)

The housing may only be opened by trained specialists!

Standard settings:


Usually the ③ IR remote control is supplied fully assembled and equipped ex works. This means that no further changes need to be made to the hardware at this point. If there are several (up to 4) bed controls in one room, it must be ensured that each remote control controls its own bed. For this, the addresses of the respective ③ IR remote controls and **BIF8plus** must be set:

- **③ Change the adress of the remote control**
Press  and  with the appropriate adress (red number).
- **Change BIF8plus adress with the ③ IR remote control**
Press  and  again and confirm with .

Address/button 1 is used as the default address in the factory settings.

Individual settings with BIF software:

- **① Connect BIF8plus to power supply**
Supply the **BIF8plus** with power (220V).
 - **Connect BIF8plus to PC**
Connect the **BIF8plus** to the PC via USB.
 - **Start software**
Then start the Universal Configurator software.
 - **Assign functions**
Now you can assign the functions to the appropriate relays by clicking on them. Confirm your final entry with **① Set**.
-

Error diagnosis activation & switching contacts:

BIF8plus recognizes when one of the switching contacts is no longer electrically separated due to a defect. In this case, the activation of the switching contacts is deactivated. For the safety of the user, further switching commands are ignored. This faulty condition is indicated by the **red LED** on the front of the housing.

10. Maintenance / Review

The device has to be checked **every 3 months**. The check must be done by a **trained specialists**. The following must be checked in detail:

- Functional test
- Connector and cable for damage and contamination
- Mounting parts for secure holding function and damage

11. Support

If you need any support or for technical questions about the product, please contact the reseller where you purchased this product. If necessary, e. g. the reseller is not available or cannot help, send us an email to support@csslabs.de. Please be sure to include a description of the problem and the serial number of the device.

12. Technical specifications

Dimensions: 145 x 106 x 45 mm

Weight: ca. 663 gr

IP Rating: IP65 (protected against water jets from all directions), RAL9002, front lugs for mounting

Power supply: Power pack 100-240 V AC, stab. 12VDC /700mA;
meets IEC 60601-1

Switch contacts: 48 V DC / 500 mA; 11 functions, 1 release contact

Infrared coding: 38 KHz, RC5 code

Temperature ranges

Use: 0 °C to 45 °C

Transportation and storage: -10 °C to 50 °C

13. Conformity/Standards

 This product is CE compliant.

14. Environmental protection and disposal

All electrical and electronic devices must not be disposed with general household waste, but must be disposed at designated disposal points. This product is subject to the European Directive 2002/96/EC and must therefore be properly disposed at the WEEE collection points. Please dispose used batteries and rechargeable batteries in a recycling bin for used batteries or via special resellers.

WEEE DE94635635

ASSISTIVE TOOLS FOR A BETTER LIFE

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© 2025 CSS MicroSystems GmbH