DE/EN

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**INSTRUCTIONS
FOR USE**

AssistX AS one

C.S.S.
MICRO SYSTEMS

Gebrauchsanweisung V3.2 DE
Instructions for Use V3.2 EN
ab / from Version V1.06
Art #20446 AssistX AS one,
Art #20076 Connect Set

CSS MicroSystems GmbH
In den Hofäckern 16
74189 Weinsberg
Germany
www.csslabs.de
© 2021 CSS MicroSystems GmbH

Inhalt

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WEEE DE94635635

Table of Contents

1. Introduction	25
2. Intended use	25
3. Explanation of symbols	26
4. Safety instructions	27
5. Cleaning and disinfection	30
6. Scope of delivery	31
7. Device Buttons, Ports, Display	32
8. Use of the device	36
9. Settings	36
Set Impulse Mode & Tone Mode:	36
Tone Mode:	37
Impulse Mode:	37
Teach in sound preset:	38
Settings of Sensitivity, tone range & hold time:	39
Discard/change settings:	39
10. Maintenance / Review	40
11. Support	41

12. Technical Data	41
13. Conformity/Standards	42
14. Environmental protection and disposal	42

1. Introduction

Thank you for choosing a high quality CSS product. Please read these instructions for use carefully and take care of these instructions when using the device. Please keep the instructions for use at hand when you use the device.

2. Intended use

AssistX AS one is a sensor that may be triggered non-contacting, only by voice. Through its varied setting options the acoustic sensor may be well adapted to the users vocal capabilities and environmental noise is effectively suppressed. **AssistX AS one** comes with an **AssistX Port** and thus may be used together with all compatible AssistX devices like **AssistX Call**, **AssistX Mobil** and **Clickman** (new generation with red sensor socket).

- Responses to noise that matches the preset volume, duration and frequency
- User friendly multi color LED for setup and feedback
- Special impuls mode for users with limited vocal capabilities
- AS one is self-monitored; In case of malfunction it triggers a call

3. Explanation of symbols

Safety instructions!

Important instructions!

4. Safety instructions

Handling

Please handle the device carefully. The device contains sensitive electronic components. It is not allowed to do any modification. The device is not life-sustaining or life-supporting. The device consists of several components. Small items may be detached from the device, which may result in suffocation or other danger to children. Small children or people with intellectual disabilities should not reach or use the device without supervision.

Do not use the device in damp, wet, very cold or very hot environments (see temperature information in technical data).

Never replace the built-in batteries by yourself! The battery can be damaged, which can lead to overheating and can lead to injuries. The integrated rechargeable lithium polymer battery may only be replaced by authorized personal.

Mounting

The device must be mounted in a way that it cannot fall down. Only recommended mounting components from the manufacturer are to be used. The user bears the sole risk of mounting.

The manufacturer is not liable for personal injury or property damage caused by the device or its components falling down.

Repair

Under no circumstances you should open the device and repair it by yourself. Disassembling the device can permanently damage the device and cause injury. Contact your reseller or the manufacturer if the device is damaged or does not work properly. If the device or its accessories are mechanically damaged, they may no longer be used.

Interaction with medical devices

The device and its components can emit electromagnetic radiation and have integrated magnets. The device must therefore not be placed on a person or brought closer than 15 cm to people using other active medical devices. The device must also not be placed directly on active medical devices. Please pay special attention to people with implanted pacemakers. If necessary, please contact the device manufacturer of the other active medical device.

Attention! Magnets

The mounting plate contains a strong neodym magnet. Magnets can affect cardiac pacemakers and implanted defibrillators.

- The pacemaker could be switched into testmode and cause discomfort
- A defibrillator could stop working
- Keep enough distance to magnets if you use such devices
- Advice other users of such device against getting too close to magnets
- Keep magnets away from all devices and objects, which can be damaged by strong magnetic fields

Recommended safety distance is 6 inch (15 cm).

Magnets cause a far-reaching magnetic field. They could damage TV sets, laptops, PC harddrives, creditcards, ATM cards, data storage, mechanical watches, hearing aids and speakers.

Vehicle keys, USB sticks, cameras, cell phones and smart phones contain non-magnetic storage media. Therefore, static magnetic fields near those devices cannot delete data. There are countless covers with magnetic closures or magnetic holders on the market, which also speaks against the notion that magnets can damage these devices. However, it cannot be ruled out that very strong magnetic fields might magnetise and maybe damage mechanical parts or built-in speakers. When in doubt, keep these devices away from strong magnets.

5. Cleaning and disinfection

Before cleaning the device, switch it off and disconnect all cables. The device should be cleaned if required with a disinfectant and cleanser intended for use in hospitals. Cleaning should only be done with a wipe. Wipe all surfaces and let them dry or follow the wiping instructions.

Do not use sprays that may allow moisture to enter the device. Under no circumstances may any component of the device be subjected to steam or wet sterilization.

6. Scope of delivery


- **AssistX AS one**
- 2 m 4 pin connection cable
- Instructions for Use
- in addition (Conection Set): connect box and power supply

WITH CONNECT SET:

7. Device Buttons, Ports, Display

- 1 LED & confirmation lamp
 - 2 AssistX Port
 - 3 Microphone
 - 4 LED display description
 - 5 Hold time (in Tone Mode)
 - 6 Sensitivity
 - 7 Tone range (in Tone Mode)
 - 8 USB programming connection
-

1 LED: show settings

The 1 LED shows when the sensor triggers or records. It has also the function of a reassurance lamp. In addition the respectively settings are shown with the associated colors. To see the current setting of sensitivity, range or hold time, press the corresponding button. When pressing a button, the 1 LED color represents the actual setting of the respective parameter. After release, the 1 LED keeps its color for 1 sec. Then it turns off.

5 Hold time (in Tone Mode)

Device in tone mode: Time span over which the sound must be kept in the set range.

	blue:	0,5 sec
	green:	1,0 sec
	red:	1,5 sec
	yellow:	2,0 sec
	white:	2,5 sec

6 Sensitivity

If the sensitivity is very high, the sensor reacts to very quiet sound.

	blue:	high
	green:	
	red:	
	yellow:	
	white:	low

7 Tone range 🎵 (in Tone Mode)

Tolerated deviation from a taught tone up and down.

 blue: 0,5 sec

 green: 1,0 sec

 red: 1,5 sec

 yellow: 2,0 sec

 white: 2,5 sec

AS one reacts on Mode:

Impulse Mode to all noise impulses fitting to the volume

Tone Mode specific to the taught in tone frequencies

8. Use of the device

Connect the device with the enclosed appropriate cable to the **AssistX Port** of your CSS product. Place the AS one sensor near the user, if necessary with an appropriate CSS holder. Because the sensor is magnetic, it adheres to all iron and steel surfaces. Follow now the steps of 9. Settings.

9. Settings

AS one can learn the frequency that triggers a sensor reaction. The sensor reacts to noise that comes with defined frequency and volume and that is held for an adjustable duration. Optional, in impulse mode the sensor reacts to any noise that reaches an adjustable volume.

Set Impulse Mode & Tone Mode:

Plug the jack cable into the **AssistX AS one**. Then press and hold the **7** tone range button for tone mode or the **5** time range button for impulse mode and then insert the plug on the other side of the cable on the **AssistX** device or **Clickman**. Press and hold the desired button for 1 sec after inserting it. As soon as the button is released, the **1** LED lights up **red** or **green** and the selected mode is active.

1 LED display

 green: Tone Mode (default)

 red: Impulse Mode

Tone Mode:

Default settings preset the device to a frequency of 300Hz +/- 100 Hz and a hold time of 1 sec. This setting covers the human voice for the most part. Sensitivity, holdtime and tolerated deviation from the preset frequency, can be adjusted using the three pushbuttons. Best detection is given for a sounding „Mmmmm“.

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Teach in sound preset:

AssistX AS one can be adapted to the users vocal capabilities. The user hums for 1.5 sec. The frequency of the sound sample is saved as the frequency basic setting. Additionally the tolerated deviation from the preset frequency can be adjusted. (tone range)

- **start recording (teach default settings)**

Press and hold the **5** time button [⌚] and **7** range button [🎵] at the same time. Release both buttons (the order is arbitrary). The device confirms with a beep.

- **recording sound**

When the sound stops, the **1** LED lights up **blue** for 1,5 sec. AS one now samples the sound that comes with its integrated microphone: *Start your humming noise now.*

- **finishing recording**

If the soundsample came with enough volume, the **1** LED lights up **green** and the frequency of the soundsample is saved. In case of a very weak soundsample the **1** LED lights up **yellow** and the Frequency presetting remains unchanged.

Settings of Sensitivity, tone range & hold time:

The settings of the AS one can be displayed and set by using the buttons. When setting the sensitivity, tone range and holding time, the ❶ LED indicates the selected setting range with a color code.

Discard/change settings:

To see the current setting of sensitivity, tone range or hold time, press the corresponding button.

- ***starting changes***

By repeated pressing of the corresponding button, you can move through the whole adjustment range. The ❶ LED shows the actual setting. 1 sec after the button is released, the ❶ LED starts blinking for 2 sec.

- ***enter changes***

Press the button while the ❶ LED is blinking to save the new setting. AS one confirms with a longer sound. The new value has been accepted.

- ***discard changes***

Leave the button released until the blinking ⓘ LED turns off after 2 sec. AS one keeps on using the last saved setting and confirms with the buzzer sounding short.

10. Maintenance / Review

The device usually does not require any maintenance. However, we recommend the following from time to time:

- functional test
 - Connectors for damage and soiling
 - Damage to the housing
-

11. Support

If you need any support or for technical questions about the product, please contact the reseller where you purchased this product. If necessary, e. g. the reseller is not available or cannot help, send us an email to support@csslabs.de. Please be sure to include a description of the problem and the serial number of the device.

12. Technical Data

Dimensions: 50 x 50 x 27 mm

Weight: 60 g (sensor without cable)

IP Rating: IP 54

Power supply: 12 VDC via **AssistX Port**

Connection: **AssistX Port**, 3.5 mm jack

(Connect Set: via power supply and connect box)

Range of frequency: 80 - 600 Hz (in Tone Mode)

Temperature ranges:

Use: 0 °C to 45 °C

Transportation and storage: -10 °C to 50 °C

13. Conformity/Standards

 This product is CE compliant.

AssistX AS one meets the requirements of VDE 0834 largely in conjunction with **AssistX Call** and **AssistX Mobil**.

14. Environmental protection and disposal

All electrical and electronic devices must not be disposed with general household waste, but must be disposed at designated disposal points. This product is subject to the European Directive 2002/96/EC and must therefore be properly disposed at the WEEE collection points. werden.

WEEE DE94635635

THE PASSIONATE ASSISTIVE TECHNOLOGY MANUFACTURER