

SWITCHES AND ACCESSORIES

Accessories

EN

Acoustic Switch (Human Voice)



AssistX AS one -
Connect Set

With Acoustic switch AS one +
AssistX Connect.

#20076

Proximity Switches



Complete Set
ProX one

Complete Set ProX one with
mounting armProximity Switch
12 mm with a mounting arm and
a universal.

#20430



Complete Set
ProX one
on tripod

Complete Set ProX one on tripod
with weight Proximity Switch
12 mm on a tripod with weight for
heavy duty.

#20428

Hygienic Pneumatic Switches



Hy-Pneumax
hinged bracket
complete set

Including pneumatic switch, fle-
xible tube, hygiene filter.

#20037



Hy-Pneumax
tripod
complete set

Including pneumatic switch, fle-
xible tube, hygiene filter.

#20035



Hy-Pneumax
6-Months
-Hygienic Pack

10x Filter, 16x mouthpiece, 16x
hose 2 m

#20031

5 Pack
hygienic filters5 pack hygiene filter suitable for
Hy-Pneumax hygiene filter with
13 cm hose and nipple.

#20036

Hy-Pneumax
mouthpiece1 Pack
5 Pack#20038
#20038-5P

packed hygienically

Hy-Pneumax
Silikonschlauch
2 m

packed hygienically

#20039

Single Buttons








KAJO Button
KB25 ø 25 mm

● red, cable 20 cm	#19131
● rot, cable 150 cm	#19132
● black, cable 20 cm	#19133
● black, cable 150 cm	#19134
○ white, cable 20 cm	#19141
○ white, cable 150 cm	#19142
● orange, cable 20 cm	#19143
● orange, cable 150 cm	#19144
● blue, cable 20 cm	#19145
● blue, cable 150 cm	#19146
● green, cable 20 cm	#19147
● green, cable 150 cm	#19148








KAJO Button
KB50 ø 50 mm

● red, cable 20 cm	#19135
● red, cable 150 cm	#19136
● black, cable 20 cm	#19137
● black, cable 150 cm	#19138
○ white, cable 20 cm	#19151
○ white, cable 150 cm	#19152
● orange, cable 20 cm	#19153
● orange, cable 150 cm	#19154
● blue, cable 20 cm	#19155
● blue, cable 150 cm	#19156
● green, cable 20 cm	#19157
● green, cable 150 cm	#19158


 PikoButton
 ø 30 mm

	blue	#19024
	yellow	#19025
	green	#19023
	orange	#19022
	red	#19020
	black	#19026
	white	#19021


 PikoButton
 ø 50 mm

	blue	#19011
	yellow	#19012
	green	#19010
	orange	#19009
	red	#19007
	black	#19015
	white	#19008







 PikoButton NC
 ø 50 mm

	red	#19007-NC
---	-----	-----------


 Big
 Buddy Button
 ø 115 mm

	red	#19100
---	-----	--------


 Buddy Button
 ø 63 mm

	blue	#19089
	yellow	#19088
	green	#19087
	red	#19086
	black	#19091

Special Switches


 KAJO Softbutton
 - blue

Soft Button blue (ø 85 mm)	#19140
Replacement Fabric Cover blue	#20351
Replacement Foam Inlay	#20213



Finger Button
FB30

Single switch

#20435



Mini Joystick

Connector D-Sub9 - for
SwitchMouse

#19092



Micro-Light-Taster

#19017

Interfaces



ClickMan

ClickMan is an programmable
adaptive sensor adapter and mul-
ti-function switch. It has an AssistX
port and a built-in battery.

#20422



Key Switch Box

With the functions Space + Enter,
2x jack socket, USB.

#20124



Mouse Switch Box
USB

With the functions left + right
mouse button, 2x jack socket,
USB.

#20127



Multi2 Switch Box

Rotary switch for 10 functions,
2 jack sockets, USB.

#20126



SwitchMouse V2
USB

Mouse simulator with USB
connection.

#20128



Double Check

Switch Tester.

#20436

Mounting Sets



Set Mounting Arm

Set Mounting Arm for AS one,
IR one, Rx one, ABL. Includes
magnetic mounting plate, Base
Clamp.

#20097



Set Mounting Arm
Kajo Button 50

Set of articulated mounting arm
Kajo Button 50, with base clamp,
also for PikoButton 50/Buddy
Button.

#20185



Set Mounting Arm
Kajo Button 25

Set of articulated mounting arm
Kajo Button 25, with base clamp,
also for PikoButton 30.

#20362



Set Mounting Arm
AssistX Mobil

AssistX Mobil articulated moun-
ting set with Rehadapt device
holder with base clamp.

#20367



Set flexible mount
AssistX Mobil

Set of flexible holder AssistX Mo-
bil, Rehadapt L3D clamp/joint for
pipes 13-35 mm.

#20368



Set Tripod
Kajo Button 50

Floor stand set Kajo Button 50,
with short gooseneck for PikoBut-
ton 50/Buddy Button.

#20369



Set Tripod
Kajo Button 25

Set of floor stand Kajo Button 25,
with short gooseneck for PikoBut-
ton 30.

#20370



Set Tripod
AssistX Mobil

Set of floor stand AssistX Mobil
with Rehadapt device holder,
ball joint.

#20371

Mounting Components



Mounting Plate
Kajo 50

Mounting Plate Kajo 50 3/8" inner
thread and 1/4" external threa-
dalso for Piko Button 50/Budy
Button.

#20360



Mounting Plate
Kajo 25

Mounting Plate Kajo 25 3/8" inner
thread and 1/4" external threadal-
so for Piko Button 30.

#20361



Gooseneck
20 cm

Gooseneck 20 cm, 11 mm diame-
ter black 3/8" and 3/8" external/
inner thread.

#19050



Base Clamp

Base Clamp for Monting Arm.

#20181



Mounting Arm

Mounting Arm, 1/4"-3/8" thread, wheel lock.

#19004



Universal Tripod

Universal Tripod with 2,5 kg stabilization weight 210/7 3/8".

#20059



Ball Joint black

Ball Joint black 1/4" external thread, incl. adapter 1/4" to 3/8".

#19032


 Tube holder
Size 2

Tube holder Size 2 for Hy-Pneumax-Module. (Exchange part)

#20432


 Tube holder
Size 1

Tube holder Size 1 for ClickMan. (Exchange part)

#20424


 Wall Mount
Size 1

Wall Mount Size 1 for ClickMan.

#19013


 Wall Mount
Size 2









Wall Mount Size 2 for Hy-Pneumax-Module.

#19014


 Mounting Plate
magnetic

Mounting Plate for AS one, IR one, Rx one, ABL. Magnetic mounting plate for Mounting Arm or wall mount.

#20459

	Wall mount magnetic	Wall mount magnetic for Remo-Box, AS one, ABL.	#20458
	Wall mount/ Mounting plate magnetic for AssistX Call	Wall mount/Mounting Plate magnetic for AssistX Call, white/ redwall mount, incl. Screws and Velcro.	#20467
	Wall mount/ Mounting plate magnetic for AssistX Repeater	Wall mount/Mounting Plate magnetic for AssistX Repeater, blackwall mount, incl. Screws and Velcro.	#20372
	Wall Mount for AssistX Mobil	Wall mount for AssistX Mobil, incl. Velcro and screws.	#20474
	Magnetic Mount for AssistX Mobil	Magnetic Mount for AssistX Mobil.	#20488
	Clamp for AssistX Mobil	Clamp for AssistX Mobil for Tube of Tripod and Mounting Arm.	#20111
	Clamp for AssistX Mobil with Base Clamp	Flexible device holder for AssistX Mobil for pipes and surfaces.	#20524
	Adapter External Thread	Thread adapter external thread 1/4" to 3/8".	#20363


 Adapter
Internal Thread

Adapter Internal Thread 1/4" to 3/8", also for Base Clamp or 3/82 external thread.

#20364


 Clamp
AssistX Mobil
Rehadapt for
Mounting Arm

Device holder for AssistX Mobil on the Mounting Arm.

#20365


 Clamp
AssistX Mobil
Rehadapt L3D
10mm

Device holder for AssistX Mobil for Rehadapt L3D 10 mm.

#20366

Drinking Aids


 Drinking Aid
complete set

Drinking Aid complete set tripod includes tripod, hose, flexible arm.

#20013


 Drinking Aid
complete set

Drinking Aid complete set flexible arm, clamp, hose, flexible arm.

#20012


 Bottle holder
for tripod

Bottle holder for tripod up to Ø 30 mm.

#19084

TABLE OF CONTENTS

Remo Environmental Control	29
Infrared Phones	29
Infrared Bed Control and IR Control Boxes	29
Easywave Radio Switching Components	31
Power Switch Adapter Sets	32
AssistX Call Systems	33
AssistX Call - Interface for Call Systems	33
AssistX Mobil - Call System	33
Interface cables for Pulseoximeters, Ventilators etc.	34
AssistX-Sensors	36
Hygienic Pneumatic Switches	36
AssistX Switches (NC)	37
Mounting AssistX	38
AssistX Accessories	40
AssistX Demo-Kits	41
AMAneo Mouse Adapter	42
Switches and Accessories	43
Acoustic Switch (Human Voice)	43
Proximity Switches	43
Hygienic Pneumatic Switches	43
Single Buttons	44
Special Switches	45
Interfaces	46
Mounting Sets	47
Mounting Components	48
Drinking Aids	51

CSS MicroSystems GmbH

In den Hofäckern 16

74189 Weinsberg

Germany

www.csslabs.de

office@csslabs.de

© 2021 CSS MicroSystems GmbH

Valid from 10/2021

For commercial customers only!

Our terms & conditions apply.

www.cssslabs.de/agb