

AssistX-CALL SYSTEM

SECURITY CONNECTS PEOPLE

#16462 V1.0 DE



Scope of delivery

The **AssistX** call system set with **transmitter**, **receiver**, and charging stations.



- **Secure calling by pressing a button:**

The **AssistX** Kajo Button with a 50 mm or 100 mm diameter.



- **Secure calling automatically when leaving bed:**

The nevisCura out-of-bed sensors – suitable for all care beds.



For care beds by **Stieglmeyer/Burmeier**, integrated out-of-bed sensors are available.

SECURE CALLING!

At Home & in Residential Groups



AssistX

Assistive Call Systems

The secure, flexible, wireless call system for home and residential groups.



CSS MicroSystems GmbH
In den Hofäckern 16
74189 Weinsberg
Germany
www.csslabs.de

© CSS 2025



assistx.de




AssistX

Assistive Call Systems

Secure calling of care givers and relatives

In home care, family members and caregivers face many challenges. Often, the person in need of care is alone in the room while the caregiver is occupied with houseworks or other tasks. The secure, wireless **AssistX call system** offers **two essential ways to trigger a call in home care situations**.



By pressing a button (1) or automatically when leaving bed: (2)
The call is sent through the transmitter of the AssistX Mobil call system.

Also for seniors and
people with dementia

AssistX-Call System - EVERY CALL IN SAFE HANDS!



The **AssistX call system** ensures that the person being cared for can **call** for help at any time – and that the caregiver or a family member securely **receives** the call, **no matter where they are in the house!**

