

AMAneo




Stops Shaking!
The Assisive Mouse Adapter.

AMAneo stops shaking!

The assistive mouse adapter AMAneo helps people with hand tremors to use the mouse on the computer or tablet much easier or to make it possible for the first time. AMAneo electronically filters the tremors of hand movements. The mouse pointer moves smoothly on the display again, even if the tremor is very strong. This is especially helpful for people with Essential Tremor and Parkinson's.

Just plug it in between the mouse and PC, notebook or iPad, and mouse control will be much easier. The intensity of the tremor filter can be adjusted. In addition, a click delay (the mouse cursor stops while mouse click) and an automatic click function are integrated.



Without tremor filter

With tremor filter

Works with almost every standard or assistive mouse. We recommend wireless mice (not Bluetooth®).



#20522 AMAneo BTi
for iPhone and iPad

Made for
iPhone | iPad



#20520 AMAneo USB
for PC, Notebook, Mac

PC/MAC/
Notebook

MADE IN GERMANY



THE PASSIONATE ASSISTIVE TECHNOLOGY MANUFACTURER

www.csslabs.de